

Windows 7

Your robot is equipped with an onboard computer. **Windows 7 Professional SP1** has been installed on this computer as its operating system (OS). This is a 32-bit OS.

All drivers and software necessary to operate your robot and its accessory devices have also been installed and configured. This includes but is not limited to:

- ARIA, core robot software development kit, with documentation, examples and source code included. (Open ARIA under the MobileRobots category in the Start menu, or C:\Program Files\MobileRobots\Aria)
- Other MobileRobots SDKs and software, depending on options ordered with the robot. This may include ARNL, SONARNL or MOGS for intelligent localization/positioning and autonomous navigation (Open ARNL, SONARNL or MOGS under the MobileRobots category in the Start menu, or C:\Program Files\MobileRobots\Arnl), and software for use with optional accessory devices ordered with the robot. (ARNL and MOGS are provided only if optional navigation packages were purchased.)
- The MobileEyes, Mapper3 and MobileSim tools. (Available in the MobileRobots category in the Start menu)
- All drivers for the onboard computer such as networking, wireless networking (wifi), video, sound, etc.

Software downloads, and information about drivers and software, as well as product manuals and a knowledge base of information about using the robot and software, are also available at <http://robots.mobilerobots.com>.

Initial Windows Setup:

The first time you boot the computer, Windows will perform some initial setup steps. Connect a keyboard, mouse and monitor to the robot the first time you turn on the computer in order to perform this initial setup. During initial setup, you will choose language and locale, enter the license key, choose a name for the system, and create a new user account. You can also perform Windows Activation at this time – see below for information on connecting the computer to your network for internet access.

The license key is printed on a sticker attached to the Windows 7 disk provided with the robot.

We recommend running Windows Update to apply all recent updates to Windows 7. You will need to connect to the internet to do so (see below). Run Windows Update by opening the Start Menu, All Programs, Windows Update.

Account Passwords:

We must create a user account when you first setup Windows. You may use this account to log after connecting keyboard, mouse and monitor, or via Remote Desktop. The Administrator account is normally disabled, but can be enabled using the Windows Administrative Tools:

1. Open the Start menu, and select **Control Panel**
2. Either open the “View By” drop-down menu, select “Small Icons” or “Large Icons”, or search for “Administrative Tools” in the search box.
3. Open **Administrative Tools**
4. Open **Local Security Policy**